


Anti-Wif1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y5W5

Item Number: CM0028753



Introduction

Reliable rabbit polyclonal antibody against human WIF1 (Wnt inhibitory factor 1) for WB, IF/ICC, ELISA; cross-reacts with mouse and rat; ideal for Wnt signaling studies.

[Learn More](#)

Parameter	Value
Product Name	WIF1 Rabbit pAb
Remarks/Alias	WIF-1; WIF1
Species	Human
GeneID (Human)	11197
GeneID	11197
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 307-377 of human WIF1 (NP_009122.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5W5
Protein Weight	42kDa
Shipping	Ice bag