

Anti-Wipf1 Polyclonal Antibody For Wb, Ihc-P, Elisa - 043516

Item Number: CM0023256



Introduction

WIPF1 Rabbit pAb validated for WB, IHC-P, ELISA. Recognizes human, mouse, rat WIPF1. Recombinant immunogen from aa 400-500. Ideal for actin cytoskeleton and immunology research.

[Learn More](#)

Parameter	Value
Product Name	WIPF1 Rabbit pAb
Remarks/Alias	WIP; WAS2; PRPL-2; WASPIP; WIPF1
Species	Human
GeneID	7456 (human)
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43516
Protein Weight	51kDa
Shipping	Ice bag