

Anti-Wipi1 Rabbit Monoclonal Antibody For Wb, Elisa - Q5Mnz9

Item Number: CM0008274

Introduction

High-specificity rabbit monoclonal antibody against human WIPI1, validated for Western blot and ELISA. Detects a 49 kDa band and cross-reacts with mouse and rat. Ideal for autophagy and cancer research.

[Learn More](#)

Parameter	Value
Product Name	WIPI1 Rabbit mAb
Remarks/Alias	ATG18; ATG18A; WIPI49; WIPI1
Species	Human
Gene ID (Human)	55062
Gene ID	55062
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5MNZ9
Protein Weight	49kDa
Shipping	Ice bag