

Anti-Wipi1 Rabbit Polyclonal Antibody For Wb, Elisa - Q5Mnz9

Item Number: CM0028704

Introduction

The WIPI1 Rabbit Polyclonal Antibody (Q5MNZ9) is validated for WB and ELISA, with cross-reactivity to human, mouse, and rat. Ideal for autophagy research and IVD raw material sourcing.

[Learn More](#)

Parameter	Value
Product Name	WIPI1 Rabbit pAb
Remarks/Alias	ATG18; ATG18A; WIPI49; WIPI1
Species	Human
Gene ID (Human)	55062
Gene ID	55062
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 187-446 of human WIPI1 (NP_060453.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5MNZ9
Protein Weight	49kDa
Shipping	Ice bag