

Anti-Wipi2 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Y4P8

Item Number: CM0027456



Introduction

Rabbit polyclonal antibody against human WIPI2, validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for studies on autophagy machinery and disease research.

[Learn More](#)

Parameter	Value
Product Name	WIPI2 Rabbit pAb
Remarks/Alias	Atg21; ATG18B; CGI-50; IDDSSA; WIPI-2; WIPI2
Species	Human
Gene ID (Human)	26100
Gene ID	26100
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 187-436 of human WIPI2 (NP_057087.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y4P8
Protein Weight	49kDa
Shipping	Ice bag