

Anti-Wisp2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - 076076

Item Number: CM0013450

Introduction

Anti-WISP2 rabbit polyclonal antibody, validated for WB, IF/ICC, ELISA, reacts with human, mouse, rat; targets CCN5 protein (27 kDa) involved in bone turnover and cell adhesion. Ideal for osteoblast signaling studies.

[Learn More](#)

Parameter	Value
Product Name	WISP2 Rabbit pAb
Remarks/Alias	CT58; WISP2; CTGF-L
Species	Human
Gene ID (Human)	8839
Gene ID	8839
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 151-250 of human WISP2 (NP_003872.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	076076
Protein Weight	27kDa
Shipping	Ice bag