

Anti-Wls Rabbit Polyclonal Antibody For Wb, Elisa - Q5T9L3

Item Number: CM0020167



Introduction

Rabbit polyclonal antibody against human WLS (GPR177/EVI), validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for Wnt signaling, protein trafficking, and developmental research. Backed by CamelBio's quality sourcing.

[Learn More](#)

Parameter	Value
Product Name	WLS Rabbit pAb
Remarks/Alias	EVI, MRP, ZKS, GPR177, mig-14, C1orf139, WLS
Species	Human
GeneID (Human)	79971
GeneID	79971
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 43-232 of human WLS (NP_001002292.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5T9L3
Protein Weight	62kDa
Shipping	Ice bag