

Anti-Wnk1 Monoclonal Antibody For Wb - Q9H4A3

Item Number: CM0004291

Introduction

Rabbit monoclonal antibody targeting human WNK1 (Q9H4A3). Validated for Western blot and ELISA; cross-reacts with mouse. Ideal for researching WNK1 kinase signaling in blood pressure regulation and ion transport.

[Learn More](#)

Parameter	Value
Product Name	WNK1 Rabbit mAb
Remarks/Alias	KDP; PSK; p65; HSN2; HSN2; PRKWNK1; PPP1R167; WNK1
Species	Human
GeneID (Human)	65125
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 1-100 of human WNK1 (Q9H4A3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9H4A3
Protein Weight	251kDa
Shipping	Ice bag