


Anti-Wnk1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9H4A3

Item Number: CM0024580



Introduction

Rabbit polyclonal antibody against human WNK1, validated for Western blot, IF/ICC, and ELISA. WNK1 is a key kinase in the WNK-SPAK/OSR1 cascade regulating blood pressure, ion homeostasis, and cell volume, essential for renal and cardiovascular research.

[Learn More](#)

Parameter	Value
Product Name	WNK1 Rabbit pAb
Remarks/Alias	KDP, PSK, p65, HSN2, HSN2, PRKWINK1, PPP1R167, WNK1
Species	Human
GeneID (human)	65125
GeneID	65125
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human WNK1 (NP_055638.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9H4A3
Protein Weight	251kDa
Shipping	Ice bag