


Anti-Wnk1 Rabbit Polyclonal Antibody For Wb And Elisa - Q9H4A3

Item Number: CM0017583



Introduction

Rabbit polyclonal antibody targeting human WNK1 (KDP, p65, HSN2), validated for Western blot and ELISA, with cross-reactivity to mouse and rat. Suitable for studies on ion homeostasis, blood pressure regulation, and cell volume response.

[Learn More](#)

Parameter	Description
Product Name	WNK1 Rabbit pAb
Remarks/Alias	KDP; PSK; p65; HSN2; HSN2; PRKWNK1; PPP1R167; WNK1
Species	Human
Gene ID (Human)	65125
Gene ID	65125
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human WNK1 (NP_055638.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9H4A3
Protein Weight	251kDa
Shipping	Ice bag