

Anti-Wnk3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Byp7

Item Number: CM0016512



Introduction

WNK3 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA. Detects human, mouse, rat WNK3. Ideal for studying electrolyte homeostasis and cell volume regulation. From CamelBio's IVD raw material portfolio.

[Learn More](#)

Parameter	Value
Product Name	WNK3 Rabbit pAb
Remarks/Alias	PRKWINK3; WNK3
Species	Human
GeneID (Human)	65267
GeneID	65267
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 500-600 of human WNK3 (NP_065973.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BYP7
Protein Weight	198kDa
Shipping	Ice bag