

Anti-Wnt1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P04628

Item Number: CM0023178

Introduction

Rabbit polyclonal antibody against human WNT1, validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for Wnt signaling, neurodevelopment, and bone biology research.

[Learn More](#)

Parameter	Value
Product Name	WNT1 Rabbit pAb
Remarks/Alias	INT1; OI15; BMND16; WNT1
Species	Human
Gene ID (Human)	7471
Gene ID	7471
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 231-370 of human WNT1 (NP_005421.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04628
Protein Weight	41kDa
Shipping	Ice bag