

Anti-Wnt10A Polyclonal Antibody For Wb And Elisa - Q9Gzt5

Item Number: CM0015837

Introduction

Rabbit polyclonal antibody against human WNT10A, suitable for WB and ELISA. Validated for human, mouse, and rat. Ideal for research in Wnt signaling, ectodermal dysplasia, and tooth development.

[Learn More](#)

Parameter	Value
Product Name	WNT10A Rabbit pAb
Remarks/Alias	OODD; SSPS; ECTD16; STHAG4; WNT10A
Species	Human
GeneID (Human)	80326
GeneID	80326
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9GZT5
Protein Weight	46kDa
Shipping	Ice bag