

Anti-Wnt16 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Ubv4

Item Number: CM0002754

Introduction

Recombinant rabbit monoclonal antibody targeting human Wnt16 (WNT16). Validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for developmental signaling and Wnt pathway studies.

[Learn More](#)

Parameter	Value
Product Name	Wnt16 Rabbit mAb
Remarks/Alias	WNT16
Species	Human
GeneID (human)	51384
GeneID	51384
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Wnt16 (Q9UBV4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBV4
Protein Weight	41kDa
Shipping	Ice bag