

Anti-Wnt2 Rabbit Monoclonal Antibody (Wb, Elisa) - P09544

Item Number: CM0006647

Introduction

Recombinant rabbit monoclonal antibody against human WNT2, validated for Western blot and ELISA. Cross-reacts with mouse and rat. Suitable for Wnt signaling and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	WNT2 Rabbit mAb
Remarks/Alias	IRP; INT1L1; WNT2
Species	Human
Gene ID (Human)	7472
Gene ID	7472
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 231-330 of human WNT2 (NP_003382.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P09544
Protein Weight	40kDa
Shipping	Ice bag