

Anti-Wnt2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P09544

Item Number: CM0028769



Introduction

Rabbit polyclonal antibody targeting human WNT2, suitable for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for Wnt signaling, embryonic development, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	WNT2 Rabbit pAb
Remarks/Alias	IRP; INT1L1; WNT2
Species	Human
Gene ID	7472
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-360 of human WNT2 (NP_003382.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09544
Protein Weight	40kDa
Shipping	Ice bag