

Anti-Wnt2B Monoclonal Antibody For Wb, Elisa - Q93097

Item Number: CM0003376



Introduction

Recombinant rabbit monoclonal antibody against WNT2B (WNT13), validated for WB and ELISA. Detects human, mouse, rat WNT2B. Ideal for Wnt signaling and lung development research.

[Learn More](#)

Parameter	Value
Product Name	WNT2B Rabbit mAb
Remarks/Alias	WNT13; WNT2B
Species	Human
GeneID (human)	7482
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-300 of human WNT2B (Q93097)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q93097
Protein Weight	44kDa
Shipping	Ice bag