

Anti-Wnt3A Monoclonal Antibody For Wb, Ip, If/Icc, If-P, Ihc-P, Elisa - P56704

Item Number: CM0009692



Introduction

Rabbit monoclonal antibody against human WNT3A protein, validated for WB, IP, IF/ICC, IF-P, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for Wnt signaling and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	WNT3a Rabbit mAb
Remarks/Alias	WNT3A
Species	Human
GeneID (Human)	89780
GeneID	89780
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P56704
Protein Weight	39kDa
Shipping	Ice bag