

Anti-Wnt3A Polyclonal Antibody For Ihc-P, Elisa - P56704

Item Number: CM0017324



Introduction

Rabbit polyclonal antibody targeting human WNT3A, validated for IHC-P and ELISA. Cross-reacts with rat. Ideal for Wnt signaling and developmental biology research. UniProt P56704.

[Learn More](#)

Parameter	Value
Product Name	WNT3A Rabbit pAb
Remarks/Alias	WNT3A
Species	Human
GeneID (human)	89780
GeneID	89780
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-352 of human WNT3A (NP_149122.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P56704
Protein Weight	39kDa
Shipping	Ice bag