

Anti-Wnt4 Polyclonal Antibody For Wb, Elisa - P56705

Item Number: CM0027524

Introduction

High-quality rabbit polyclonal antibody against human WNT4, validated for WB and ELISA. Reacts with human, mouse, and rat. Ideal for studying urogenital and lung development.

[Learn More](#)

Parameter	Value
Product Name	WNT4 Rabbit pAb
Remarks/Alias	WNT-4; SERKAL; WNT4
Species	Human
GeneID (Human)	54361
GeneID	54361
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P56705
Protein Weight	39kDa
Shipping	Ice bag