

Anti-Wnt5A Monoclonal Antibody For Wb, Ihc-P, Elisa - P41221

Item Number: CM0002578

Introduction

Recombinant rabbit monoclonal antibody against human WNT5A. Validated for WB, IHC-P, ELISA. Specific for human WNT5A, molecular weight 42kDa. Suitable for Wnt signaling and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	WNT5A Rabbit mAb
Remarks/Alias	hWNT5A, WNT5A
Species	Human
GeneID	7474
Immunogen	Synthetic peptide (aa 281-380 of human WNT5A, P41221)
Source	Rabbit
Category	Monoclonal Antibody
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P41221
Protein Weight	42kDa
Shipping	Ice bag