

Anti-Wnt5A Rabbit Polyclonal Antibody For Wb - P41221

Item Number: CM0017412

Introduction

Rabbit polyclonal antibody targeting human WNT5A (aa 116-215), validated for WB, IHC-P, IF/ICC, ELISA. Recognizes human, mouse, rat WNT5A, a key Wnt signaling ligand involved in development and cancer.

[Learn More](#)

Parameter	Value
Product Name	WNT5A Rabbit pAb
Remarks/Alias	hWNT5A; WNT5A
Species	Human
GeneID (human)	7474
GeneID	7474
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 116-215 of human WNT5A (NP_003383.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P41221
Protein Weight	42kDa
Shipping	Ice bag