

Anti-Wnt5B Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9H1J7

Item Number: CM0028869

Introduction

Rabbit polyclonal antibody targeting human WNT5B (SWISS: Q9H1J7). Validated for WB, IHC-P, ELISA in human, mouse, rat. Ideal for Wnt signaling and developmental studies.

[Learn More](#)

Parameter	Value
Product Name	WNT5B Rabbit pAb
Remarks/Alias	WNT5B
Species	Human
GeneID (Human)	81029
GeneID	81029
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human WNT5B (NP_110402.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H1J7
Protein Weight	40kDa
Shipping	Ice bag