


Anti-Wnt7A Rabbit Polyclonal Antibody For Wb And Elisa - O00755

Item Number: CM0025796



Introduction

Rabbit polyclonal antibody against human WNT7A, validated for WB and ELISA with cross-reactivity to mouse. Targets Wnt-7a, a secreted ligand in canonical Wnt/beta-catenin signaling, crucial for embryonic development and CNS angiogenesis.

[Learn More](#)

Parameter	Value
Product Name	WNT7A Rabbit pAb
Remarks/Alias	Wnt-7a; WNT7A
Species	Human
GeneID (human)	7476
GeneID	7476
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 101-200 of human WNT7A (NP_004616.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O00755
Protein Weight	39kDa
Shipping	Ice bag