

Anti-Wnt9A Rabbit Polyclonal Antibody For Wb And Elisa - O14904

Item Number: CM0025473

Introduction

CamelBio's anti-WNT9A rabbit polyclonal antibody, validated for WB and ELISA, recognizes human, mouse, and rat WNT9A, a Wnt signaling ligand critical for bone development and joint integrity.

[Learn More](#)

Parameter	Value
Product Name	WNT9A Rabbit pAb
Remarks/Alias	WNT14; WNT9A
Species	Human
Gene ID (Human)	7483
Gene ID	7483
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 250-350 of human WNT9A (NP_003386.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14904
Protein Weight	40kDa
Shipping	Ice bag