


## Anti-Wrap53 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Bur4

Item Number: CM0020653



### Introduction

Rabbit polyclonal antibody against human WRAP53 (TCAB1), for WB and ELISA, cross-reacts with mouse. Ideal for telomere biology, Cajal body, and DNA repair research.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | WRAP53 Rabbit pAb                                                                                                  |
| Remarks/Alias    | DKCB3; TCAB1; WDR79; WRAP53                                                                                        |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 55135                                                                                                              |
| GeneID           | 55135                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 62-160 of human WRAP53 (NP_060551.2) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, ELISA                                                                                                          |
| Cross Reactivity | Mouse                                                                                                              |
| SWISS            | Q9BUR4                                                                                                             |
| Protein Weight   | 59kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |