

Anti-Wrn Rabbit Polyclonal Antibody For Wb, Elisa - Q14191

Item Number: CM0025599



Introduction

Rabbit polyclonal antibody against human WRN, validated for WB and ELISA. Targets the C-terminal region (aa 1223-1432) of the Werner syndrome helicase/exonuclease, ideal for DNA repair and genome stability research.

[Learn More](#)

Parameter	Value
Product Name	WRN Rabbit pAb
Remarks/Alias	RECQ3; RECQL2; RECQL3; WRN
Species	Human
GeneID (human)	7486
GeneID	7486
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1223-1432 of human WRN (NP_000544.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q14191
Protein Weight	162kDa
Shipping	Ice bag