

Anti-Wsb1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y6I7

Item Number: CM0013988

Introduction

Rabbit polyclonal antibody against human WSB1, validated for WB and ELISA. Recognizes WD repeat and SOCS box-containing protein 1, a component of E3 ubiquitin ligase complex involved in proteasomal degradation. - Q9Y6I7

[Learn More](#)

Parameter	Value
Product Name	WSB1 Rabbit pAb
Remarks/Alias	SWIP1; WSB-1; WSB1
Species	Human
GeneID (human)	26118
GeneID	26118
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 182-421 of human WSB1 (NP_056441.6)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9Y6I7
Protein Weight	47kDa
Shipping	Ice bag