

Anti-Wstf Rabbit Monoclonal Antibody For Wb And Elisa - Q9Uig0

Item Number: CM0008715



Introduction

Rabbit monoclonal anti-WSTF antibody for WB and ELISA applications. Recognizes human, mouse, and rat WSTF. Suitable for chromatin remodeling and DNA damage response studies.

[Learn More](#)

Parameter	Value
Product Name	WSTF Rabbit mAb
Remarks/Alias	WSTF; WBSCR9; WBSCR10
Species	Human
GeneID (Human)	9031
GeneID	9031
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human WSTF (Q9UIG0)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UIG0
Protein Weight	171kDa
Shipping	Ice bag