


Anti-Wt1 Polyclonal Antibody For Wb, Ihc-P, Elisa - P19544

Item Number: CM0015445



Introduction

Rabbit polyclonal antibody targeting human WT1, validated in WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for studying transcription regulation and cancer biology.

[Learn More](#)

Parameter	Value
Product Name	WT1 Rabbit pAb
Remarks/Alias	GUD; AWT1; WAGR; WT-1; WT33; NPHS4; WIT-2; WT1
Species	Human
GeneID (Human)	7490
GeneID	7490
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-302 of human WT1 (NP_001185480.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19544
Protein Weight	49 kDa
Shipping	Ice bag