


Anti-Wt1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P19544

Item Number: CM0005541



Introduction

WT1 Rabbit Monoclonal Antibody (P19544) validated for WB, IHC-P, ELISA. Recognizes human, mouse, rat WT1 protein. Recombinant immunogen. Suitable for urogenital development and tumor suppressor studies.

[Learn More](#)

Parameter	Value
Product Name	WT1 Rabbit mAb
Remarks/Alias	GUD; AWT1; WAGR; WT-1; WT33; NPHS4; WIT-2; WT1
Species	Human
GeneID (Human)	7490
GeneID	7490
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 60-250 of human WT1 (P19544).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19544
Protein Weight	49kDa
Shipping	Ice bag