


Anti-Wt1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P19544

Item Number: CM0021473



Introduction

WT1 rabbit polyclonal antibody for WB, IF/ICC, and ELISA, reactive in human, mouse, and rat. Targets the Wilms tumor protein, a key transcription factor in urogenital development and tumorigenesis.

[Learn More](#)

Parameter	Value
Product Name	WT1 Rabbit pAb
Remarks/Alias	GUD; AWT1; WAGR; WT-1; WT33; NPHS4; WIT-2; WT1
Species	Human
Gene ID (Human)	7490
Gene ID	7490
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-302 of human WT1 (NP_001185480.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19544
Protein Weight	49kDa
Shipping	Ice bag