

Anti-Wtap Polyclonal Antibody For Wb, If/Icc, Elisa - Q15007

Item Number: CM0018094



Introduction

Rabbit polyclonal antibody targeting human WTAP, validated for Western blot, immunofluorescence, and ELISA. Detects human, mouse, and rat WTAP, suitable for m6A methylation and splicing regulation research.

[Learn More](#)

Parameter	Value
Product Name	WTAP Rabbit pAb
Remarks/Alias	Mum2; WTAP
Species	Human
Gene ID	9589
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 212-355 of human WTAP (NP_004897.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15007
Protein Weight	44 kDa
Shipping	Ice bag