

## Anti-Wwc1 (Kibra) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8Ix03

Item Number: CM0022480



### Introduction

High-quality rabbit polyclonal antibody targeting human WWC1 (KIBRA), validated in WB, IHC-P, and ELISA. Reacts with human, mouse, and rat samples, suitable for research on Hippo signaling and cognition.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | WWC1 Rabbit pAb                                                                                                    |
| Remarks/Alias    | KIBRA; HBEBP3; HBEBP36; MEMRYQTL; PPP1R168; WWC1                                                                   |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 23286                                                                                                              |
| GeneID           | 23286                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 290-426 of human WWC1 (NP_056053.1). |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, IHC-P, ELISA                                                                                                   |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | Q8IX03                                                                                                             |
| Protein Weight   | 125kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                            |