

Anti-Wwc2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q6Awc2

Item Number: CM0020710

Introduction

Affinity-purified rabbit polyclonal antibody targeting human WWC2 (BOMB) protein, validated for WB, IHC-P, and ELISA. Detects human, mouse, and rat WWC2, a regulator of the Hippo pathway. Shipped on ice.

[Learn More](#)

Parameter	Value
Product Name	WWC2 Rabbit pAb
Remarks/Alias	BOMB; WWC2
Species	Human
Gene ID (Human)	80014
Gene ID	80014
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 160-430 of human WWC2 (NP_079225.5)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6AWC2
Protein Weight	134kDa
Shipping	Ice bag