


Anti-Wwox Polyclonal Antibody For Wb, Elisa - Q9Nzc7

Item Number: CM0015406



Introduction

Rabbit polyclonal antibody against human WWOX, validated for Western blot (WB) and ELISA. Cross-reacts with mouse and rat. Targets the WW domain-containing oxidoreductase, a tumor suppressor involved in apoptosis and Wnt signaling.

[Learn More](#)

Parameter	Value
Product Name	WWOX Rabbit pAb
Remarks/Alias	FOR; WOX1; DEE28; EIEE28; FRA16D; SCAR12; HHCMA56; PRO0128; SDR41C1; D16S432E; WWOX
Species	Human
GeneID (Human)	51741
GeneID	51741
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NZC7
Protein Weight	47kDa
Shipping	Ice bag