


Anti-Wwox Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nzc7

Item Number: CM0016815



Introduction

Rabbit polyclonal antibody targeting human WWOX (Q9NZC7), a tumor suppressor and oxidoreductase involved in apoptosis and Wnt signaling. Validated for WB, IF/ICC, ELISA; cross-reacts with mouse. Molecular weight: 47kDa.

[Learn More](#)

Parameter	Value
Product Name	WWOX Rabbit pAb
Remarks/Alias	FOR; WOX1; DEE28; EIEE28; FRA16D; SCAR12; HHCMA56; PRO0128; SDR41C1; D16S432E; WWOX
Species	Human
GeneID (Human)	51741
GeneID	51741
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human WWOX (NP_057457.1))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NZC7
Protein Weight	47kDa
Shipping	Ice bag