

Anti-Wwp1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9H0M0

Item Number: CM0028114



Introduction

Rabbit polyclonal antibody against WWP1, a NEDD4-like E3 ubiquitin-protein ligase involved in ubiquitination and proteasomal degradation. Suitable for WB and ELISA, reacts with human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	WWP1 Rabbit pAb
Remarks/Alias	AIP5; Tiul1; hSDRP1; WWP1
Species	Human
GeneID (Human)	11059
GeneID	11059
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human WWP1 (NP_008944.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H0M0
Protein Weight	105kDa
Shipping	Ice bag