

## Anti-Wwp2 Rabbit Polyclonal Antibody For Wb And Elisa - 000308

Item Number: CM0024031



### Introduction

WWP2 Rabbit Polyclonal Antibody validated for WB and ELISA. Targets human WWP2 (NEDD4-like E3 ubiquitin-protein ligase), a key regulator of ubiquitination and proteasomal degradation in stem cells and T-cells.

[Learn More](#)

| Parameter        | Value                                                                                                                |
|------------------|----------------------------------------------------------------------------------------------------------------------|
| Product Name     | WWP2 Rabbit pAb                                                                                                      |
| Remarks/Alias    | AIP2, WWp2-like, WWP2                                                                                                |
| Species          | Human                                                                                                                |
| Gene ID          | 11060 (human)                                                                                                        |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 151-400 of human WWP2 (NP_001257383.1) |
| Source           | Rabbit                                                                                                               |
| Category         | Polyclonal Antibodies                                                                                                |
| Application      | WB, ELISA                                                                                                            |
| Cross Reactivity | Human                                                                                                                |
| SWISS            | 000308                                                                                                               |
| Protein Weight   | 99kDa                                                                                                                |
| Shipping         | Ice bag                                                                                                              |