

Anti-Xanthine Dehydrogenase/Oxidase (Xdh) Monoclonal Antibody For Wb, Ihc-P, Elisa - P47989

Item Number: CM0007100



Introduction

Rabbit monoclonal antibody targeting human Xanthine Dehydrogenase/Oxidase (XDH), validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Ideal for purine metabolism and ROS generation studies.

[Learn More](#)

Parameter	Value
Product Name	Xanthine Oxidase (XDH) Rabbit PolymAb®
Remarks/Alias	XO; XOR; XAN1
Species	Human
GeneID (human)	7498
GeneID	7498
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1000-1100 of human Xanthine Oxidase (XDH). (NP_000370.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P47989
Protein Weight	146kDa
Shipping	Ice bag