

Anti-Xanthine Oxidase (Xdh) Rabbit Monoclonal Antibody For Wb - P47989

Item Number: CM0004239



Introduction

High-quality rabbit monoclonal antibody targeting human Xanthine Oxidase (XDH), validated for WB, IHC, IF, and ELISA. Ideal for purine metabolism and reactive oxygen species studies.

[Learn More](#)

Parameter	Value
Product Name	Xanthine Oxidase (XDH) Rabbit mAb
Remarks/Alias	XO; XOR; XAN1; Xanthine Oxidase (XDH)
Species	Human
Gene ID (Human)	7498
Gene ID	7498
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 202-293 of human Xanthine Oxidase (XDH) (P47989)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P47989
Protein Weight	146kDa
Shipping	Ice bag