

Anti-Xanthine Oxidase (Xdh) Rabbit Polyclonal Antibody For Wb And Elisa - P47989

Item Number: CM0013590



Introduction

High-quality rabbit polyclonal antibody targeting human Xanthine Oxidase (XDH), validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for studying purine metabolism and oxidative stress.

[Learn More](#)

Parameter	Value
Product Name	Xanthine Oxidase (XDH) Rabbit pAb
Remarks/Alias	XO; XOR; XAN1; Xanthine Oxidase (XDH)
Species	Human
Gene ID (Human)	7498
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 985-1084 of human Xanthine Oxidase (XDH) (NP_000370.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P47989
Protein Weight	146kDa
Shipping	Ice bag