

Anti-Xbp-1S Rabbit Monoclonal Antibody For Wb, Chip, Elisa - P17861-2

Item Number: CM0009132

Introduction

Anti-XBP-1s rabbit monoclonal antibody validated for WB, ChIP, and ELISA. Cross-reacts with human and mouse XBP-1s. Ideal for ER stress, UPR, and cancer research. - P17861-2

[Learn More](#)

Parameter	Value
Product Name	XBP-1s Rabbit mAb
Remarks/Alias	XBP2; TREB5; XBP-1; TREB-5
Species	Human
Gene ID (Human)	7494
Gene ID	7494
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 188-287 of human XBP-1s (NP_001073007.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ChIP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P17861-2
Protein Weight	40 kDa
Shipping	Ice bag