

Anti-Xcp1 Rabbit Polyclonal Antibody For Wb And Elisa - B6Tic7

Item Number: CM0024196

Introduction

Rabbit polyclonal antibody against maize XCP1 (cysteine protease 5) for Western blot (WB) and ELISA applications. Recognizes maize XCP1, predicted 71 kDa, ideal for plant protease research.

[Learn More](#)

Parameter	Value
Product Name	XCP1 Rabbit pAb
Remarks/Alias	IDP105, GRMZM2G066326, XCP1, Cysteine protease5
Species	Maize
GeneID	100283433
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 157-377 of maize XCP1 (NP_001149806.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Maize
SWISS	B6TIC7
Protein Weight	71kDa
Shipping	Ice bag