


Anti-Xcr1 Polyclonal Antibody For Ihc-P And Elisa - P46094

Item Number: CM0020330



Introduction

Rabbit anti-human XCR1 polyclonal antibody, validated for IHC-P and ELISA, with rat cross-reactivity. Targets chemokine XC receptor 1 (lymphotactin receptor) involved in chemokine signaling. Suitable for immunodetection applications.

[Learn More](#)

Parameter	Value
Product Name	XCR1 Rabbit pAb
Remarks/Alias	GPR5; CCXCR1; XCR1
Species	Human
GeneID (Human)	2829
GeneID	2829
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human XCR1 (NP_001019815.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Rat
SWISS	P46094
Protein Weight	39kDa
Shipping	Ice bag