


Anti-Xiap Rabbit Monoclonal Antibody (Ko Validated) For Wb, Elisa - P98170

Item Number: CM0002253



Introduction

KO-validated rabbit monoclonal antibody targeting human XIAP (BIRC4/IAP3), with cross-reactivity to mouse and rat. Ideal for Western blot and ELISA to study apoptosis, NF-kB signaling, and ubiquitination pathways.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] XIAP Rabbit mAb
Remarks/Alias	API3; ILP1; MIHA; XLP2; BIRC4; IAP-3; hIAP3; hIAP-3; AP
Species	Human
GeneID (Human)	331
GeneID	331
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human XIAP (NP_001158.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P98170
Protein Weight	57kDa
Shipping	Ice bag