


Anti-Xiap Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P98170

Item Number: CM0018707



Introduction

Rabbit polyclonal antibody targeting human XIAP (X-linked inhibitor of apoptosis protein). Validated for WB, IF/ICC, ELISA in human, mouse, rat. MW ~57kDa. For apoptosis, NF-κB signaling, and ubiquitination studies.

[Learn More](#)

Parameter	Value
Product Name	XIAP Rabbit pAb
Remarks/Alias	API3; ILP1; MIHA; XLP2; BIRC4; IAP-3; hIAP3; hIAP-3; XIAP
Species	Human
GeneID (Human)	331
GeneID	331
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human XIAP (NP_001158.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P98170
Protein Weight	57kDa
Shipping	Ice bag