

Anti-Xk Polyclonal Antibody For Wb, If/Icc, Elisa - P51811

Item Number: CM0017966

Introduction

Rabbit polyclonal antibody targeting XK (XKR1, Kx antigen). Suitable for Western blot, immunofluorescence, and ELISA in human, mouse, rat. Investigate ER membrane adapter protein function and VPS13A recruitment.

[Learn More](#)

Parameter	Value
Product Name	XK Rabbit pAb
Remarks/Alias	KX; NA; NAC; X1k; XKR1; XK
Species	Human
Gene ID (Human)	7504
Gene ID	7504
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 59-138 of human XK (NP_066569.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51811
Protein Weight	51kDa
Shipping	Ice bag