


Anti-Xpc Rabbit Polyclonal Antibody For Wb - Q01831

Item Number: CM0024964



Introduction

Rabbit polyclonal antibody targeting human XPC (XP3/RAD4/p125) protein. Validated for Western blot, IHC-P, IP, and ELISA. Ideal for nucleotide excision repair and transcription coactivation studies. SWISS: Q01831.

[Learn More](#)

Parameter	Value
Product Name	XPC Rabbit pAb
Remarks/Alias	XP3; RAD4; XPCC; p125; XPC
Species	Human
Gene ID	7508
Immunogen	Recombinant protein containing amino acids 700-940 of human XPC (NP_004619.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human
SWISS	Q01831
Protein Weight	106kDa
Shipping	Ice bag