


Anti-Xpd/Ercc2 Monoclonal Antibody For Wb, Ihc-P, Elisa - P18074

Item Number: CM0002103



Introduction

Anti-XPD/ERCC2 rabbit monoclonal antibody for Western Blot, immunohistochemistry (paraffin), and ELISA. Cross-reacts with human, mouse, rat. Ideal for studying transcription-coupled nucleotide excision repair and TFIIH complex function. Predicted molecular weight 87kDa.

[Learn More](#)

Parameter	Value
Product Name	XPD/ERCC2 Rabbit mAb
Remarks/Alias	EM9; TTD; XPD; TTD1; COFS2; TFIIH; XPD/ERCC2
Species	Human
GeneID (Human)	2068
GeneID	2068
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 600-700 of human XPD/ERCC2 (P18074)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18074
Protein Weight	87kDa
Shipping	Ice bag